

REGO-LOK™

For Securing CGA Fittings on Liquid Cylinders

The REGO-LOK™ is designed for installation on the Rego T9450 and T9460 Series liquid cylinder valves to deter and prevent the removal of the CGA fitting from the valve. The REGO-LOK™ retains standard CGA outlet connection fittings so unauthorized persons do not remove the fitting. By use of a special one-way bolt, the REGO-LOK™ is secured to the valve.

The REGO-LOK™ installs in a few minutes with the use of screwdrivers. It installs without valve disassembly, brazing, welding, or drilling. The REGO-LOK™ deters and prevents fitting removal by gas customers, however allows the replacement of fillings by authorized gas supplier plant personnel.

Use REGO-LOK™ for compliance with CGA SB-26 for medical and industrial liquid cylinders.

Satisfies CGA SB-26 and FDA requirements for medical and industrial liquid cylinders.

Stainless Steel REGO-LOK™ with one-way bolt.

Retrofit all common liquid cylinder valves.

Can be supplied on new Rego liquid cylinder valves.

REGO-LOK™ indicates *WARNING: DO NOT REMOVE*

Worn CGA fittings can be simply replaced by authorized personnel (Requires new 9464RL-6 Bolt)

Can fit over existing fittings for CGA 540, CGA 440, CGA 295, CGA 320, and CGA 326. Check fitting hex size.

Prevents loosening of CGA fittings on valves.

NOTE: Rego supplied fitting PIN CGA5B0RL is required for REGO-LOK™ use with CGA 580 connection.

Part Number	Item Description	Typical Service Connection
9464RL-0	For 3/4" hex fittings	CGA295
9464RL-1	For 1/4" fittings	CGA440, CGA320 & CGA326
9464RL-2	For 1" fittings	CGA540
9464RL-3	For 1 1/4" hex CGA5B0RL fitting by Rego	CGA560
CGASB0RL	3/4" MNPT x CGA for use with 9464RL-3	CGA580

REGO-LOK™ is a registered trademark of RegO Patent applied for.